



VISION SAFETY FUNNELS

VISIONEC-LAB | WWW.VISIONEC.HU

Up to 5 capillary connection (1/16", 2mm, 3mm, 1/8") from HPLC for continuous solvent disposal



Filtered air, pressure release



Liquid



Standard 1/4"-28 ports

If the funnel is used without HPLC, the port can be sealed with blind plug. Exhaust filter must be used to avoid pressure build up in the container

Manual valve made of PTFE to provide exceptional chemical resistance

Harmful vapors



Liquid from funnel and from capillaries

Wetted materials are:
PP (Polipropylene)
PTFE (Teflon)

Safety Funnel with shut off valve for DIN/S61 containers



Troughput: 30ml/sec
(water@25°C)

Dimensions

Cap size: DIN61
Total height: 185mm
Trough hole dimension: 6x12mm
Funnel volume: 350ml
Funnel width: 74mm@top

Part No.

Description

VLF61-0002-TW

S60/S61 Safety Cap with safety funnel and 2 ports;
1 x 1/4"-28 fitting for 1/8" tubing + open port for exhaust filter

VLF61-0011-TW

S60/S61 Safety Cap with safety funnel and 2 ports;
1 x barbed connector for 5-8mm ID tubing + open port for exhaust filter

Safety Funnel with shut off valve for DIN45/GL45 containers



Troughput: 25ml/sec
(water@25°C)

Dimensions

Cap size: GL45/DIN45
Total height: 140mm
Trough hole dimension: 5x8mm
Funnel volume: 100ml
Funnel width: 48mm@top

Part No.

Description

VLF45-0002-TW

DIN45 Safety Cap with safety funnel with 2 ports;
1 x 1/4"-28 fitting for 1/8" tubing + open port for exhaust filter

VLF45-0011-TW

DIN45 Safety Cap with safety funnel with 2 ports;
1 x barbed connector for 5-8mm ID tubing + open port for exhaust filter



VisionEC products are available worldwide. For the distributors in your country, contact VisionEC-Lab Hungary, Corporate Office at info@visionec.hu



Terms and Conditions

Subject to VisionEC-Lab Standard Terms & Conditions, which may be viewed at www.visionec.hu

© 2018 VisionEC-Lab Ltd. All rights reserved.